

Hgb Pro

Professional Hemoglobin Testing System®

CLIAwaived.com™



For more information please contact CLIAwaived Inc. at
11578 Sorrento Valley Road ■ San Diego, CA 92121

Tel: (858) 481-5031 ■ Toll Free: (888) 862-7739

E-mail: info@cliawaived.com

Visit us on the web at: www.cliawaived.com

Distributed by: CLIAwaived, Inc.



Precision and accuracy can be affordable.

Trust the point-of-care industry leaders.

SIMPLE TO USE

The Hgb Pro Professional Hemoglobin Testing System is portable and easy to use. Simply insert a test strip into the meter and place a small drop of blood onto the strip. Since the Hgb Pro is self-calibrating, no additional calibration is required.

FEATURES

- 250 memory data storage
- Large, easy-to-read LCD screen
- Operates on one 3-volt lithium battery
- Cost-effective

- CLIA-waived
- Small blood-drop sample
- Results within 30 seconds
- Lab-quality accuracy and precision

SYSTEM OVERVIEW

THREE SIMPLE STEPS

Step 1: Turn on meter.

Step 2: Insert test strip.

Step 3: Apply blood.

QUALITY CONTROL

Liquid quality controls, Para 4[®] (product number 215110), are available from Streck[®] Laboratories (1-800-843-0912 or www.streck.com) for the Hgb Pro. It is recommended that quality control testing be performed upon receipt of a new lot of test strips and once every three months thereafter.

ORDERING INFORMATION

HGBPROKIT: Starter kit (Contains one meter and one bottle of 25 test strips.)

HGBMETER: One Hgb Pro meter

HGBTEST100: Test strips (Contains four bottles of 25 test strips.)

REIMBURSEMENT INFORMATION

Use CPT code 85018QW.

Hgb Pro Professional Hemoglobin Testing System[®] is a registered trademark of ITC. Streck[®] and Para 4[®] are registered trademarks of Streck Laboratories.



8 Olsen Avenue • Edison, NJ 08820 USA
tel: 732.548.5700 • fax: 732.248.1928
www.itcmed.com

a subsidiary of Thoratec Corporation

